



IN DCS- CORDANT HEALTH SOLUTIONS TREATMENT PROVIDER TRAINING

©2020 Cordant Health Solutions®. All rights reserved. Do not
distribute without written permission from Cordant Health Solutions.

Cordant
Health Solutions®

AGENDA

- Changes in DCS Process
- Who is Cordant Health Solutions?
- The Collection Process
 - Oral Fluid Collection
 - Monitored Urine Collections
- Sentry
- Reading Results
- Additional Resources

IN DCS SPECIFIC CHANGES IN PROCESS

What is Changing?

- Changing from current vendor to Cordant Health Solutions
- Urine specimen collection should be “monitored” rather than “observed”
 - If an adulteration event is found or other concerns warranted- full observation will be required for that individual
- Provider will have the option to utilize either oral fluid or urine at their discretion
- Provider can request additional panels if warranted

What is not Changing?

- SUDA and SUDT referrals will still come with the units of testing that you can administer for assessments and treatment
- Treatment providers will have access to the results in the Cordant System and the results will also be sent by Cordant directly to KidTraks

CORDANT HEALTH SOLUTIONS

Beyond our commitment to excellence and providing the best solutions available, we are driven by a higher purpose at Cordant Health Solutions—to make a difference in the lives of those we serve, especially those affected by addiction. **We are all part of the solution.**

Cordant's National Footprint

DEDICATION TO SERVICE


- **4.7 million samples** tested annually across 50 states
- **8,000 patients served** by our managed pharmacy program
- **65,000 MAT prescriptions** filled
- **395,000 pills disposed** of through our take-back program

LEADER IN QUALITY

- **1 of only 21 labs** nationwide certified by the U.S. Department of Health and Human Services (SAMHSA)
- CLIA, CAP and CAP FDT certified


INDUSTRIES WE SERVE

- Substance Use Treatment
 - *Abstinence-Based, Medication-Assisted*
- Criminal Justice Agencies
 - *Treatment Courts, Parole, Probation, Social Services*
- Mental Health
- Hospitals and Health Systems
- Chronic Pain Management
- Workers' Compensation



Laboratory Services Available Nationally

Pharmacy Services Available Regionally



THE COLLECTION PROCESS

BEFORE YOU COLLECT



1. Verify the name and date of birth of the client or patient
2. Search and Print collection order from Sentry
3. Review Process Acknowledgment form while explaining the collection process
4. If any issues occur during the collection, “report violation” or notes on section page
5. Samples are to always remain in clear sight to both the client and the collector until the chain of custody is complete and then sealed and packaged to send to the lab



ORAL FLUID COLLECTION




ORAL FLUID SPECIMEN COLLECTION

COLLECTING THE SPECIMEN

1. Always ensure the client has not used tobacco products, drank any fluids or consumed any food products within the last **10 minutes** before collection. This includes gum, candies and mouthwash.
2. Always inspect the expiration date of the device and open package- discard if expired
3. Have the client grab the device handle and position the cotton pad under their tongue
4. The client must not speak or bite the device during the collection process
5. The collector must remain in the room with the client until the collection is complete
6. The collection device must stay under the tongue until the Volume Adequacy Indicator turns **BLUE**
7. The indicator should turn blue in 2-5 minutes but could take up to 10
 - ❖ **What happens if the indicator does not turn blue after 10 minutes?**
 - ❖ Discard the device
 - ❖ Instruct the client to drink water
 - ❖ Then wait 10 minutes
 - ❖ Recollect with new device

If indicator does not turn blue, you will not have enough oral fluid to ensure accurate testing- Resulting in a potential FALSE NEGATIVE




ORAL FLUID SPECIMEN COLLECTION

SECURING THE SPECIMEN

1. Have the client firmly hold the transport tube upright
2. Remove the red cap by pushing up ward with the thumb
3. Be sure none of the liquid spills or splashes out of tube
4. Have the client grab the handle of the collection device and remove from their mouth
5. Place the device bad end first into the transport tube
6. Replace the red cap and push down until it SNAPS into place

❖ If device pad is bitten off, falls off or is in some way damaged, discard sample and device and recollect

- Placing the device in upside down will be rejected by lab!
- Potential false negative because drug won't be properly extracted from the pad into the solution





MONITORED URINE COLLECTIONS

MONITORED URINE SPECIMEN COLLECTION

COLLECTING AND SECURING THE SPECIMEN

1. Blue septic agent should be utilized before every monitored collection event
2. Ask the client to remove coats/jackets, purses, etc. Have them remove all pocket contents; personal items are to remain outside of the collection restroom
3. Be sure to check around the client's waist band and body looking for any baggies, strings, devices, etc that seem out of place. Have the client wash their hands thoroughly (no soap) to thwart any attempts to conceal items on or around the hands or arms
4. Instruct client not run any water or flush until specimen is handed back to collector
5. **Monitored collections:** The collector stands outside of the door to listen for any potential improprieties (Collector is NOT observing the client urinating into the container)
6. **Observed collections:** Inform the client you must observe the urine leaving the body
7. Lab requires a minimum of 15 mL of specimen
8. Visually check sample for signs of contamination or tampering (note any visual observations in the Test request space provided on COC)
9. Check temperature on temp strip. This should be between 90 and 100 degrees. Record on form. If temperature did not register, recollect and notate first failed attempt in Sentry

Cordant

CONFIDENTIAL

SENTRY CHAIN OF CUSTODY

Cordant

Forensic Solutions™

Test Request & Chain of Custody Document

P.O Box 70,000
Flagstaff, AZ 86003

AL1092190

REQUESTED BY: Nearpass, Ben (10712) PRINTED BY: Jorgensen (Norahem), Brandon (8692) COLLECTED AT: Technical Resource Management	
TRM	
TRM, INC. NEW HIRE P.O. BOX 70000 800-348-4422 FLAGSTAFF, AZ 86004	
I certify that the specimen accompanying this form is mine, and was sealed with a tamper evident seal in my presence and that the information provided on this form and on the label is correct.	
_____ DONOR SIGNATURE	01/09/2020 DATE
I certify that the specimen accompanying this form is the specimen given to me by the donor. I further certify that it was sealed in my presence with the accompanying specimen evidence seal.	
_____ COLLECTOR SIGNATURE	01/09/2020 DATE
(ONLY COMPLETE THIS SECTION IF NEEDED) _____ ADDITIONAL CUSTODIAN/OTHER THAN COLLECTOR	
DATE	

DONOR NAME (OR SPECIMEN I.D.) LAST FIRST M.I.				
SSN		CLIENT ID # 123456789		OTHER ID 8158
DONOR TELEPHONE DAY 6024787954				
DATE COLLECTED 01/09/2020	TIME COLLECTED 13:29 MST	DATE OF BIRTH n/a	SEX ☒ M ☐ F	RACE
LIST PRESCRIPTION & OVER THE COUNTER MEDICATIONS TAKEN IN THE LAST 10 DAYS 			TEMP (90-100 F) ☐	
			VISUALLY MONITORED ☐	
Specimen inspected for: *SEAL INTEGRITY *MATCHING COC # *DONOR INITIALS				
TEST REQUEST: 280 - 8 Drug Screen "includes Etg"			PLACED IN TEMPORARY STORAGE	
			COC defects listed here: (if none are listed, none were found)	

SEND TO LAB

AL1092190

AL1092190

COC #:AL1092190
Date Collected:01/09/2020
Time Collected:13:29 MST

Donor Name
Donor Case #: 123456789
Donor SSN:
Donor ID: 8158

Tests Requested
280 - 8 DRUG SCREEN "INCLUDES ETG"

DONOR RECEIPT

DONOR NAME (OR SPECIMEN I.D.)				
LAST		FIRST		M.I.
SSN		CLIENT ID #	OTHER ID	
		123456789	8158	
DONOR TELEPHONE DAY 6024787954				
DATE COLLECTED	TIME COLLECTED	DATE OF BIRTH	SEX	RACE
01/09/2020	13:29 MST	n/a	☒ M ☐ F	
LIST PRESCRIPTION & OVER THE COUNTER MEDICATIONS TAKEN IN THE LAST 10 DAYS			TEMP (90-100 F) ☐	
			VISUALLY MONITORED ☐	
			Specimen inspected for: *SEAL INTEGRITY *MATCHING COC # *DONOR INITIALS	

Donor demographics will be automatically populated within Sentry-
No manual entry necessary

Cordant
Forensic Solutions™

Document

P.O Box 70,000
Flagstaff, AZ 86003

AL1092190

REQUESTED BY: Nearpass, Ben (10712) PRINTED BY: Jorgensen (Norohem), Brandon (8692) COLLECTED AT: Technical Resource Management	DONOR NAME (OR SPECIMEN I.D.) LAST FIRST M.I.
TRM  TRM, INC. NEW HIRE P.O. BOX 70000 800-348-4422 FLAGSTAFF, AZ 86004	SSN CLIENT ID # 123456789 OTHER ID 8158
DONOR DAY 6024787954	

I certify that the specimen accompanying this form is mine and was sealed with a tamper evident seal in my presence and that the information provided on this form and on the label is correct.			
DONOR SIGNATURE	01/09/2020 DATE		
I certify that the specimen accompanying this form is the specimen given to me by the donor. I further certify that it was sealed in my presence with the accompanying specimen evidence seal.			
COLLECTOR SIGNATURE	01/09/2020 DATE		
(ONLY COMPLETE THIS SECTION IF NEEDED)			
ADDITIONAL COLLECTOR OTHER THAN COLLECTOR	DATE		

DATE COLLECTED 01/09/2020	TIME COLLECTED 13:29 MST	DATE OF BIRTH n/a	SEX ☒ M ☐ F
LIST PRESCRIPTION & OVER THE COUNTER MEDICATIONS TAKEN IN THE LAST 10 DAYS 			TEMP (90-100 F) ☐
TEST REQUEST: 280 - 8 Drug Screen "Includes Etc"			VISUALLY MONITORED ☐ Specimen Inspected for: *SEAL INTEGRITY *MATCHING COC # *DONOR INITIALS
			PLACED IN TEMPORARY STORAGE
COC defects listed here: (If none are listed, none were found)			
By: _____			Date: _____

SEND TO LAB

SEND TO LAB

AL1092190



AL1092190

COC #:AL1082190
Date Collected:01/09/2020
Time Collected:13:29 MST

Donor Name: _____
Donor Case #: 123456789
Donor SSN: _____
Donor ID: 8158

Tests Requested
280 - 8 DRUG SCREEN "INCLUDES ETG"

DONOR RECEIPT

01/09/2020

DONOR SIGNATURE

DATE

01/09/2020

COLLECTOR SIGNATURE

DATE _____

- Donor and Collector signature
- Tamper evident seal

PLACE SEAL
ACROSS CAP
OF SPECIMEN
CONTAINER

Donor
Initials _____

Date _____




PLACE
OVER
CENTER
OF CAP


E2362073

CORDANT REQUISITION FORM

The diagram shows the layout of the Cordant Requisition Form. It includes a section for 'Donor Initials' and 'Date' with lines for handwritten input. A central circular area is labeled 'PLACE OVER CENTER OF CAP'. To the right is a barcode with the number 'E1963235' printed below it.

Once label is adhered, have client initial and date the label




Ensure seal is placed correctly over the sample lid and down the sides of the vial or device

- This acknowledges the sample was sealed in the client's presence and is indeed their sample
- Seal number **MUST** match the number on the requisition form

PROPER CHAIN OF CUSTODY

FINAL CHECK: Both the collector and the client must sign the requisition form



Make sure the bar codes, initials, name and signature all match!

PROCESS ACKNOWLEDGMENT FORM

- Completed with each collection
- Sent to the lab with COC and sample
- Acknowledges the collection protocol was followed and Chain of Custody was maintained
- Urine, oral fluid and hair version

Chain of Custody
ment

1 BLACK BALLPOINT PEN ONLY / NO HIGHLIGHTERS!

3 NAME (OR SPECIMEN ID) PRINT LEGIBLY

CASE # OTHER ID #

NOR DAY ()

E2362073

Cordant
Health Solutions

Process Acknowledgment Form – Oral Fluid Sample

Indiana DCS has requested that you provide a sample for drug testing. Upon collection, your sample will be sealed and labeled with a chain of custody sticker and sent to Cordant for testing. The results will be provided to Indiana DCS. The following collection process will be followed:

1. The collector will visually inspect your mouth to ensure it is empty.
2. It is expected you wait 10 minutes after eating, drinking, smoking, chewing gum, etc. before the oral fluid collection.
3. The collector will provide you with the collection device. Place the collection swab under your tongue. Do not chew on the device, talk or drink beverages during collection.
4. Once the indicator turns blue, you will place the device swab side down into the collection vial.
5. You will initial a chain of custody label.
6. You will sign and date the chain of custody form.
7. The sample will then be sealed with the chain of custody label and placed in a tamper evident bag in your presence.

Review and sign

I acknowledge the above process was followed as outlined.

Donor Signature _____ Date _____ Time _____

Signature of Guardian if donor is a minor _____ COC # _____

Collector Signature _____ Date _____ Time _____

42 CFR part 2 prohibits unauthorized disclosure of these records.

RELEASE OF INFORMATION FORM

- Completed with each collection
- Sent to the lab with COC and sample
- This release ensures compliance with a contractual requirement between Cordant and DCS - all requirements of HIPAA and 42 CFR Part 2 must be maintained.

CONSENT TO RELEASE HEALTH AND SUBSTANCE ABUSE RECORDS

1. I, _____ (hereinafter, the "Individual"), by and through this consent to release health and substance abuse records (this "Consent and Release") hereby authorize Technical Resource Management, LLC d/b/a Cordant Health Solutions (the "Testing Entity") to receive, use, release, disclose, and re-disclose the following health and substance use information (the "Health and Substance Abuse Records"), by and among themselves, and the Indiana Department of Child Services (DCS):
- ☐ The entirety of all medical, substance abuse, and other records relating to the Individual maintained by the Testing Entity
- ☐ Other (please explain): _____
2. I authorize the Testing Entity to receive, use, release, disclose and re-disclose the above-described information by and among themselves, and to the Indiana Department of Child Services (DCS), for the purposes of securing, administering, documenting, and processing child welfare services and benefits by and through DCS and in responding to subpoenas and court orders for such records in relation to legal matters involving DCS and the administration of DCS benefits and services pending before the State and Federal Courts in the State of Indiana, unless otherwise specified below:
- ☐ Other (please specify): _____
3. The authorization to receive, use, release, disclose, and re-disclose of the above-described Health and Substance Abuse Records shall expire two (2) years after the date of signature hereto or upon the following occurrence:
- ☐ Please specify event, if not a specific date (ex., Closure of DCS Case Number:): _____
4. I understand that my Health and Substance Abuse Records are protected under federal law, including the final regulations governing the confidentiality of substance use disorder patient records, 42 C.F.R. Part 2, and the Health Insurance Portability and Accountability Act of 1996 ("HIPAA"), 45 C.F.R. Parts 160 and 164, and cannot be disclosed without my written consent pursuant to this Consent and Release unless otherwise provided within their implementing regulations. Federal law prohibits the person or organization to whom disclosure is made from making any further disclosure of this information unless further disclosure is expressly permitted by the written authorization of the person to whom it pertains or as otherwise permitted by 42 C. F. R. Part 2.
5. I understand that if the person or entity to whom the Testing Entity may disclose or re-disclose the above-described Health and Substance Abuse Records is not a doctor, health care provider or health plan, the Health and Substance Abuse Records may not be protected by HIPAA, and that recipient may use or disclose the Health and Substance Abuse Records to other non-covered entities. I understand that the Health and Substance Abuse Records may include information relating to sexually transmitted diseases, acquired immunodeficiency syndrome (AIDS), or human immunodeficiency virus (HIV); behavioral or mental health services; and treatment for alcohol and drug abuse.
6. I understand that my refusal to sign this Consent and Release will not affect the Individual's ability to obtain services from the Testing Entity; however it may affect the Individual's ability to secure or otherwise participate in child welfare services and benefits administered through DCS which may require access to the Individuals Health and Substance Abuse Records.
7. I understand I have the right to inspect or copy the Health and Substance Abuse Records disclosed pursuant to this Consent and Release and to receive a copy of the signed Consent and Release. I understand that I may revoke (cancel) this Consent and Release at any time and that such revocation must be in writing and delivered to the Testing Entity. The Testing Entity cannot be held responsible for having disclosed Health and Substance Abuse Records in reliance upon this Consent and Release before receiving a written revocation.
8. I understand that the Testing Entity and their respective workforce members are released from legal responsibility or liability for disclosing Health and Substance Abuse Records authorized by this Consent and Release.
9. I acknowledge that I understand and had an opportunity to ask questions before I signed, this Consent and Release.

Signature of Individual (or Individual's Personal Representative)

Date of Signature

Print Name of Individual or Personal Representative (if applicable)

Individual's Date of Birth

*Personal Representative Only (please describe your authority to sign on behalf of the Individual (e.g., Parent, Guardian, etc.)):

42 CFR part 2 prohibits unauthorized disclosure of these records.

SHIPPING TO LAB



SEAL IN BAG


1. Place specimen in the smaller compartment of the bag with the absorbent pad
2. Fold both the Process Acknowledgement Form, Release of Information Form and the Chain of Custody in half and then in half again
3. Place the folded forms in the larger compartment of the specimen bag, separate from the specimen
4. Fold the adhesive flap over the cross-hatch slit opening.
5. Place in Fed Ex Box or Bag for shipping to the lab




SENTRY



HIGH LEVEL PROCESS



SENTRY LOG-IN SENT THROUGH ENCRYPTED EMAIL




- For SUD Providers that have more than one location, a different log-in will be provided for each location.
- Click on the link in the email
- You will be prompted to either create a password OR enter the one you've already created.

IMPORTANT: Sentry log-In credentials will not be sent until we have received the signed Qualified Service Organization Agreement (QSOA).

LOGIN

- Login at <https://sentry.cordanthhs.com/Sentry>




Primary Functions within Sentry for SUD Providers:

- ❖ View/download results
- ❖ Perform collections

FINDING THE DCS CLIENT IN SENTRY

- There are two ways to find the DCS client:
 1. Type the name of the client in the quick search box
 2. Click on the “Granted Clients” tab. All DCS clients that you have access to will be included on the Granted Access tab.
- If you cannot find the client, please call our referral management team.
 - The referral management team will need to verify that we have actually received the referral.
 - *Important note: Cordant will only receive the drug testing referral after you have accepted the treatment referral*



CLIENT PROFILE IN SENTRY

SENTRY
Cordant Health Solutions

Quick Search

Client Info

Test Results

Testing Schedule

Phone System (IVR) Log

Collections

Client History

client abc, test

Edit

You have been granted access to this client as of 11/07/2019 by My Training Account of Chris Stephens Training.
[remove access]

Name: client abc, test

Client ID#: 7656757

User: Jr, CSTRAIN

Group: 2x / week + 10% 7 Drug

SSN:

Date of Birth: N/A

Age: N/A

Sex: M

Race: N/A

Compliance Score: 100 / 100 / 100 / 100 / 100 [?]

Last Result Date: N/A

Next Scheduled Test: N/A

Last Phone System (IVR) call-in: N/A [See call log]

Phone System (IVR) Code: 550030

Active Dates: n/a

Prescription Medications: [Edit]

Photos

Upload

Webcam Capture

No images found

Notes

Add note

Date Entered	Re: Date	Negativity	Author	Message
07/09/2020 11:07 AM		5	Stephens (Norchem), Chris	Client moved to 2x / week + 10% 7 Drug(92295)
11/08/2019 2:46 PM		5	Account, My Training	Created 11/08/2019 by Account, My Training (cstrain)

Alerts

- Once you click on the client's name, the client's profile in Sentry will pull up.
- All of the tabs that are displayed across the top of the screen relate to this specific client.
- The two primary tabs that will be utilized are the Test Results tab and the Collections tab (highlighted with the RED arrows)

THE ABILITY TO REMOVE GRANTS TO REDUCE CLUTTER

OPTIONAL FEATURE

SENTRY™
Cordant Health Solutions

Quick Search

Client Info | Test Results | Testing Schedule | Phone System (IVR) Log | Collections | Client History

client abc, test [Edit](#)

Name: client abc, test
Client ID#:: 7656757
User: Jr, CSTRAIN
Group: 2x / week + 10% 7 Drug
SSN:
Date of Birth: N/A **Age:** N/A
Sex: M **Race:** N/A

Compliance Score: 100 / 100 / 100 / 100 / 100 [?]
Last Result Date: N/A
Next Scheduled Test: N/A
Last Phone System (IVR) call-in: N/A [See call log]
Phone System (IVR) Code: 550030
Active Dates: n/a
Prescription Medications: [\[Edit\]](#)

Photos [Upload](#) [Webcam Capture](#)
No images found


Notes [Add note](#)

Date Entered	Re: Date	Negativity	Author	Message
07/09/2020 11:07 AM		5	Stephens (Norchem), Chris	Client moved to 2x / week + 10% 7 Drug(92295)
11/08/2019 2:46 PM		5	Account, My Training	Created 11/08/2019 by Account, My Training (cstrain)

- The referral received by Cordant includes the parent name of the Treatment Provider
 - There could be multiple locations for the provider
 - Access to the DCS client will be granted to all log-ins for the SUD Provider.
- An individual location can click on the “remove access” button to remove clients that are not being treated at that particular location.
 - The “remove grant” feature would only apply to the specific log-in that is removing the grant.


TEST RESULTS

- The “Test Results” tab will provide you with all results for the DCS client that you have access to.
- Click on the hyper link in the “Accession #” column to view the final result report
- You can also view the completed COC form by clicking on that hyper link.
- The data can be sorted by clicking on any of the column headers.
- **Important note: Cordant will also provide ALL TEST RESULTS to Indiana DCS.**



Test Date	Lab Received	Result Released	Accession #	COC #	Result	Creatinine	Abnormal Reason(s)
07/20/2018	07/20/2018	07/20/2018	VL0035919	AJ720239B	Issue / Abnormal - Not Forgiven		Violation Attempt
03/06/2018	03/06/2018	03/06/2018	VL0032371	AJ306129D	Issue / Abnormal - Not Forgiven		Violation Attempt
09/07/2016	09/07/2016	09/07/2016	VL0020804	AG9070E9C	Issue / Abnormal - Not Forgiven		Violation Attempt
09/06/2016	09/06/2016	09/06/2016	VL0020777	AG9061A86	Issue / Abnormal - Not Forgiven		Violation Attempt
07/27/2016	07/27/2016	07/27/2016	VL0020006	AG7270D2E	Issue / Abnormal - Forgiven		Violation Attempt
01/11/2016	01/11/2016	01/11/2016	VL0016838	AG1110C1D	Issue / Abnormal - Forgiven		Violation Attempt
11/10/2015	11/10/2015	11/10/2015	VL0015918	AFB100AA8	Issue / Abnormal - Not Forgiven		Violation Attempt
05/11/2015	05/11/2015	05/11/2015	VL0013580	AF5110660	Issue / Abnormal - Not Forgiven		Violation Attempt
04/08/2015	04/08/2015	04/08/2015	VL0013260	AF408037B	Issue / Abnormal - Not Forgiven		Violation Attempt
01/15/2016		01/15/2016	FTK97540		Abnormal - Not Forgiven		THC/Marijuana,Opiates,Methamphetamine
01/15/2016		01/15/2016	FTK97539		Normal		
08/06/2018		08/06/2018	FTK968987		Abnormal - Not Forgiven		THC/Marijuana
01/12/2016		01/12/2016	FTK95722		Abnormal - Not Forgiven		Benzodiazepines,Methamphetamine
06/28/2018		06/28/2018	FTK928181		Abnormal - Not Forgiven		THC/Marijuana,Cocaine
06/28/2018		06/28/2018	FTK928180		Abnormal - Not Forgiven		Breath Alcohol
06/03/2015		06/03/2015	FTK9246		Issue / Normal		Shy Bladder
04/09/2015		04/09/2015	FTK918		Abnormal - Not Forgiven		Amphetamine,Breath Alcohol
05/14/2018		05/14/2018	FTK881913		Abnormal - Not Forgiven		Barbiturates,Buprenorphine

DOWNLOAD TEST RESULTS FEATURE




- The download test results feature allows you to download results using a number of different formats for all OR a selection of results for the individual client.
- If you download into a PDF format, you will receive a copy of the final result report and you have the option to include the Chain of Custody in the download.

Report Options:

Filter: All Results
Abnormal only
Normal Only
Select Individual tests

Format: Microsoft Excel (.xls)
Comma-delimited (.csv)
Pipe-delimited (.psv)
Adobe Acrobat (.pdf)

TEST RESULT EMAIL NOTIFICATIONS



SENTRY™
Cordant Health Solutions

Quick Search

[Home](#) | [Settings](#) | [Log Out](#)
Technical Resource Management
Amanda Gibbs

[Your Account](#) | [Printing](#) | **[Notifications](#)**

Alerts

Notification Setup

All Results:
 E-Mail to send upon receiving any result
 E-Mail to send to nightly for all Tests

Missed IVR Call-Ins:
☐ Suppress alerts in web interface
 E-Mail to send to nightly for all callin violations

No-Show for Testing
☐ Suppress alerts in web interface
 E-Mail to send to nightly for all violations

Abnormal Tests
☐ Suppress alerts in web interface
 E-Mail to send to when result posts
 E-Mail to send to nightly for all abnormal Tests


[Update](#) | [Cancel](#)

Copyright © Cordant Health Solutions, All Rights Reserved
Phone: (800) 348-4422
Hours: Mon-Fri, 5:30am-5pm PST
clientservices@cordanths.com

[Contact Cordant](#)

Click on “Settings” and then the “Notifications” tab to set-up result notifications to be sent to a specific email address.

COLLECTIONS WILL ONLY BE ALLOWED WHEN THERE ARE VALID REFERRALS FROM INDIANA DCS



SENTRYDV™
Click to return to your home screen
Cordant Health Solutions

Quick Search

Home | Settings | Log Out
Chris Stephens Training
Customer ABC

Client Info | Test Results | Testing Schedule | Phone System (IVR) Log | **Collections** | Client History

Collecting For: client abc, test

Name: test client abc
Date of Birth: Age: N/A
Sex: M
Race: N/A
Client ID#: 7656757

There are no valid referrals

- If you see this box that there are no valid referrals, **please contact the FCM**. This means that the referral was cancelled. The FCM will need to send a new referral for us to be able to provide drug testing services for the client.
- If there are no valid referrals, the system will not allow you to proceed through the collection process.

COLLECTION PROCESS

SINGLE REFERRAL EXAMPLE

- SUD Providers will perform testing under one of two referral types:
 - Substance Use Disorder Assessment
 - Substance Use Outpatient Treatment
- Click on the “Identity Confirmed” button to proceed to the next step in the collection process.

SENTRY™
Cordant Health Solutions

Search: cynthia nix

Home | Settings | Log Out
103688 Indiana DCS
Treatment Provider

Client Info | Test Results | Testing Schedule | Phone System Log | **Collections** | Client History

Collecting For: Doe, John

Name: John Doe
Date of Birth: 04/27/1978 Age: 43
Sex: M
Race: N/A
Personal ID #: 00123456789

Selected Referral
Billable Unit: RF0006321234
Description: SUBSTANCE USE DISORDER ASSESSMENT-Oral Initial
Quantity Available: 12
Quantity Already Collected: 2

IDENTITY CONFIRMED
IDENTITY INVALID

The number of drug testing units available on the referral and the number already collected will be displayed here.

If additional drug testing units are needed, please contact the FCM.

COLLECTION PROCESS

MULTIPLE REFERRAL EXAMPLE

SENTRY™
Cordant Health Solutions

Home | Settings | Log Out
103688 Indiana DCS
Treatment Provider

Client Info | Test Results | Testing Schedule | Phone System Log | **Collections** | Client History

Collecting For: Doe, Jane

Name: Jane Doe
Date of Birth: 11/12/1995 Age: 25
Sex: F
Race: N/A
Personal ID #: 12345678910

Pick A Referral Type

☒ RF0006321234	Billable Unit: RF0006321234 Description: SUBSTANCE USE OUTPATIENT TREATMENT-Urine Initial Quantity Available: 12 Quantity Already Collected: 0
☐ RF0006332106	Billable Unit: RF0006332106 Description: RANDOM DRUG TESTING-Oral Initial Quantity Available: 26 Quantity Already Collected: 0

IDENTITY CONFIRMED
IDENTITY INVALID

If there are multiple referrals for the DCS client, it is very important to choose either the SUBSTANCE USE DISORDER ASSESSMENT referral OR the SUBSTANCE USE OUTPATIENT TREATMENT referral (as applicable) prior to selecting “identity confirmed.”

COLLECTION PROCESS

UPDATE MEDICATIONS – OPTIONAL FEATURE

- You can choose to add medications that the client is taking. These medications will print on the Chain of Custody form.
- Click “Proceed with Collection” to move to the next screen.

The screenshot displays the SENTRY web application interface. At the top, the SENTRY logo and 'Cordant Health Solutions' are visible. A search bar contains 'indiana dcs'. Navigation links for 'Home', 'Settings', and 'Log Out' are in the top right, along with 'Technical Resource Management' and 'Amanda Gibbs'. A horizontal menu includes 'Case Info', 'Test Results', 'Testing Schedule', 'IVR Log', 'Collections', and 'Case History'. The 'Collections' tab is active, showing 'Collecting For: Rushing, Charity'. Below this, a section titled 'UPDATE MEDICATIONS' contains client details: Name: Charity Rushing, Date of Birth: 09/11/1977 Age: 43, Sex: F, Race: N/A, and Case #: 2157262914. A green 'PROCEED WITH COLLECTION' button is centered. Below the button is a 'MEDICATIONS' section with a text box stating 'No medications associated with this case'. It includes a search bar 'Find medication by name...' and a date picker for 'Prescription expiration:'. An 'INCLUDE MEDICATION' button is at the bottom right of the medication section. The footer contains copyright information, contact details, and a 'Contact Cordant' link.

COLLECTION PROCESS

SELECT TESTS TO BE ORDERED

SENTRY™
Cordant Health Solutions

indiana

Case Info | Test Results | Testing Schedule | IVR Log | **Collections** | Case History

Collecting For: Rushing, Charity

UPDATE TESTS

Test Groups:

- 7200 - URINE: 9 DRUG PANEL PLUS ADULTERANTS
- 7201 - ORAL FLUID: 9 DRUG PANEL
- 1541 - URINE: KRATOM PANEL
- 3701 - URINE: SPICE PANEL
- 3601 - URINE: DESIGNER STIMULANT PANEL BY LC/MS/MS
- TAC-5870-H-CP1-E - HAIR: 10 DRUG PANEL

Hold CTRL to select multiple test groups

Addon Tests:

- 19 - URINE: PCP BY IA
- 29.1 - URINE: : ETHYL GLUCURONIDE-ETG (500 NG/ML)
- 11.7 - ORAL FLUID: : ETHANOL (0.02%)

Hold CTRL to select multiple tests

Copyright © Cordant Health Solutions, All Rights Reserved

Phone: (800) 348-4422
Hours: Mon-Fri, 5:30am-5pm PST
nclientservices@cordanthsh.com

- Tests to be ordered are selected on this screen. Hold the CTRL button to select multiple tests.
- **Two primary panels:**
 - 7200 – Urine: 9 drug panel plus adulterants
 - 7201 – Oral Fluid: 9 drug panel
- **Both of these panels include the following substances:**
 - Amphetamines/ Methamphetamines, Benzodiazepines, Buprenorphine, Cocaine, Cannabinoids, Fentanyl, Opiates, Oxycodone, and Tramadol.

AD-HOC/ADDITIONAL TEST REQUESTS

- The standard panel for both urine and oral fluid includes the following substances:
Amphetamines/ Methamphetamines, Benzodiazepines, Buprenorphine, Cocaine, Cannabinoids, Fentanyl, Opiates, Oxycodone, and Tramadol.
- Add-on substances that are included on the Cordant contract:

Urine		Oral Fluid
Alcohol (EtG)	Meperidine	Alcohol (Ethanol)
Methadone	Dextromethorphan	Methadone
PCP	Gabapentin	PCP
Barbiturates	Zolpidem	Barbiturates
Ecstasy	Spice	Ecstasy
EDDP (methadone metabolite)	Bath Salts	
6-AM (heroin metabolite)	Kratom	

If you need to test a substance that is not available for selection in Sentry, please contact Cordant so that the test selection can be added for future collections.

COLLECTION PROCESS

PRINT THE CHAIN OF CUSTODY FORM

Alerts

SENTRY™
Cordant Health Solutions

indiana dcs

Case Info

Test Results

Testing Schedule

IVR Log

Collections

Case History

COC #AM6222379

06/22/21 4:03 PM

Gibbs, Amanda (amandag)

PDF Image

COC not displayed? [click here](#) to download it.

Cancel Collection

Report Violation

Done

Report

1 / 1

89%

Cordant
Forensic Solutions™

Test Request & Chain of Custody

Document

P.O. Box 70,000
Flagstaff, AZ 86003


AM6222379

REQUESTED BY: Master, IN DCS Region 11 (41800)
PRINTED BY: Gibbs, Amanda (7007)
COLLECTED AT: Technical Resource Management

TRM



TRM, INC. NEW HIRE
P.O. BOX 70000
800-348-4422
FLAGSTAFF, AZ 86004

I certify that the specimen accompanying this form is mine,
and was sealed with a tamper evident seal in my presence

DONOR NAME (OR SPECIMEN I.D.)

LAST

FIRST

M.I.

Rushing

Charity

SSN

PERSONAL ID #

OTHER ID

2157262914

551974

DONOR TELEPHONE

DAY

(765) 557-4446

DATE COLLECTED

TIME COLLECTED

DATE OF BIRTH

SEX

RACE

Copyright © Cordant Health Solutions, All Rights Reserved

Phone: (800) 348-4422
Hours: Mon-Fri, 5:30am-5pm PST
clientservices@cordanthhs.com

Print the chain of custody and test request form on the specific paper that was sent to you. Place paper in tray and click print. The barcode should land on the security seal.

Cordant

CONFIDENTIAL

36

RESULT INTERPRETATION

URINE DRUG TESTING FINAL REPORT

Testing Results				
Test	Result	Outcome	Cutoff	Notes
Screening Tests by IA				
URINE: Meth/Amphetamines		negative	1000 ng/ml	
URINE: Benzodiazepines		negative	300 ng/ml	
URINE: Cocaine		negative	300 ng/ml	
URINE: Methadone		SEE BELOW	300 ng/ml	
URINE: Opiates		SEE BELOW	300 ng/ml	
URINE: THC		negative	50 ng/ml	
URINE: Oxycodone		negative	300 ng/ml	
URINE: Ethyl Glucuronide-ETG		SEE BELOW	500 ng/ml	
URINE: Buprenorphine/Suboxone		negative	5 ng/ml	
URINE: Fentanyl		negative	2.0 ng/ml	
URINE: Soma/Carisoprodol		negative	100 ng/ml	
URINE: Gabapentin		negative	1000 ng/ml	
Confirmation Tests				
Methadone by LC/MS/MS				
URINE: Methadone metabolite	1020 ng/ml	**POSITIVE	300 ng/ml	
URINE: Methadone	1375 ng/ml	**POSITIVE	300 ng/ml	
Opiates by LC/MS/MS				
URINE: Codeine	Not Detected	negative	300 ng/ml	
URINE: Morphine	Not Detected	negative	300 ng/ml	
URINE: Hydrocodone	380 ng/ml	**POSITIVE	300 ng/ml	
URINE: Hydromorphone	250 ng/ml	Detected	300 ng/ml	• Results consistent with Hydrocodone use. • Detected (drug fingerprint confirmed) above the laboratory limit of quantitation.
URINE: Oxycodone	Not Detected	negative	300 ng/ml	
URINE: Oxymorphone	Not Detected	negative	300 ng/ml	
ETG/ETS by LC/MS/MS				
URINE: Ethyl Glucuronide	3354 ng/ml	**POSITIVE	500 ng/ml	
URINE: Ethyl Sulfate	2184 ng/ml	**POSITIVE	100 ng/ml	• ETG AND ETS INTERPRETIVE GUIDELINES: A) Negative for ETG and ETS means negative for alcohol consumption B) POSITIVE levels for ETG and ETS are consistent with alcohol consumption C) ETG negative or Unable to Confirm and ETS POSITIVE are consistent with alcohol consumption D) ETS negative, Not Detected and ETG POSITIVE may not indicate alcohol consumption.
Specimen Validity Tests				
Specimen Validity Panel				
URINE: Creatinine	47.0 mg/dl		20 mg/dl	
URINE: Basic Adulteration Check		normal		• Specimen checked for unusual color, physical characteristics and abnormal instrument response.



Specimen Information	
Donor Name:	Collected: 04/09/2019 00:00:00
Donor DOB:	Received: 04/12/2019 12:06:57
Accession:	Reported: 04/15/2019 19:20:06
COC:	Donor Other ID:
Type (Matrix): Urine	Donor Case:
Client Code:	
Client:	
Requested By:	

Test	Result	Outcome	Cutoff	Notes
Screening Tests by IA				
URINE: Meth/Amphetamines		negative	1000 ng/ml	
URINE: Benzodiazepines		SEE BELOW	300 ng/ml	
URINE: Cocaine		negative	300 ng/ml	
URINE: Methadone		negative	300 ng/ml	
URINE: Opiates		negative	300 ng/ml	
URINE: THC		negative	50 ng/ml	
URINE: Oxycodone		negative	300 ng/ml	
URINE: Ethyl Glucuronide-ETG		negative	500 ng/ml	• ETG negative screen suggests alcohol was not consumed.
URINE: Buprenorphine/Suboxone		negative	5 ng/ml	
URINE: Fentanyl		negative	2.0 ng/ml	
URINE: Soma/Carisoprodol		negative	100 ng/ml	
URINE: Gabapentin		negative	1000 ng/ml	
Confirmation Tests				
Benzodiazepines by LC/MS/MS				
URINE: Oxazepam	Not Detected	negative	100 ng/ml	
URINE: Nordiazepam	Not Detected	negative	100 ng/ml	
URINE: Temazepam	Not Detected	negative	100 ng/ml	
URINE: Lorazepam	Not Detected	negative	100 ng/ml	
URINE: Alprazolam-Metabolite	232 ng/ml	**POSITIVE	100 ng/ml	• Results consistent with Xanax (Alprazolam) use.

All important donor demographics entered on the chain of custody will be located here to the top left corner



Specimen Information

Donor Name:	Collected: 04/09/2019 00:00:00
Donor DOB:	Received: 04/12/2019 12:06:57
Accession:	Reported: 04/15/2019 19:20:06
COC:	Donor Other ID:
Type (Matrix): Urine	Donor Case:
Client Code:	
Client:	
Requested By:	

Testing Results

Test	Result	Outcome	Cutoff	Notes
Screening Tests by IA				
URINE: Meth/Amphetamines		negative	1000 ng/ml	
URINE: Benzodiazepines		SEE BELOW	300 ng/ml	
URINE: Cocaine		negative	300 ng/ml	
URINE: Methadone		negative	300 ng/ml	
URINE: Opiates		negative	300 ng/ml	
URINE: THC		negative	50 ng/ml	
URINE: Oxycodone		negative	300 ng/ml	
URINE: Ethyl Glucuronide-ETG		negative	500 ng/ml	• ETG negative screen suggests alcohol was not consumed.
URINE: Buprenorphine/Suboxone		negative	5 ng/ml	
URINE: Fentanyl		negative	2.0 ng/ml	
URINE: Soma/Carisoprodol		negative	100 ng/ml	
URINE: Gabapentin		negative	1000 ng/ml	
Confirmation Tests				
Benzodiazepines by LC/MS/MS				
URINE: Oxazepam	Not Detected	negative	100 ng/ml	
URINE: Nordiazepam	Not Detected	negative	100 ng/ml	
URINE: Temazepam	Not Detected	negative	100 ng/ml	
URINE: Lorazepam	Not Detected	negative	100 ng/ml	
URINE: Alprazolam-Metabolite	232 ng/ml	**POSITIVE	100 ng/ml	• Results consistent with Xanax (Alprazolam) use.

- All screening information at top of result and grouped together at the drug class level
- Confirmation data, if performed, will be located at the bottom of the result
 - Each drug within drug class will be broken out
- Helpful additional interpretive comments located in the right column

Confirmation Tests				
Benzodiazepines by LC/MS/MS				
URINE: Oxazepam	Not Detected	negative	100 ng/ml	
URINE: Nordiazepam	Not Detected	negative	100 ng/ml	
URINE: Temazepam	Not Detected	negative	100 ng/ml	
URINE: Lorazepam	Not Detected	negative	100 ng/ml	
URINE: Alprazolam-Metabolite	232 ng/ml	**POSITIVE	100 ng/ml	• Results consistent with Xanax (Alprazolam) use.
URINE: Flurazepam-Metabolite	Not Detected	negative	100 ng/ml	
URINE: Clonazepam-Metabolite	Not Detected	negative	100 ng/ml	
URINE: Midazolam	Not Detected	negative	100 ng/ml	
URINE: Triazolam-Metabolite	Not Detected	negative	100 ng/ml	
Specimen Validity Tests				
Specimen Validity Panel				
URINE: Creatinine	26.5 mg/dl		20 mg/dl	
URINE: Basic Adulteration Check		normal		• Specimen checked for unusual color, physical characteristics and abnormal instrument response.
Additional Comments				
• Testing performed at Cordant Forensic Solutions, 1760 E Route 66, P.O. Box 70000, Flagstaff, AZ 86004. • Tests performed under CAP-FDT certification. • Specimen received sealed and intact unless otherwise noted. • CLIA #03D0936918, CAP-FDT #6913001 • Report Released By: KJM - B.S., Certifying Scientist.				

Validity testing for urine located on the very bottom of the report

CHAIN OF CUSTODY ERROR NOTES

Requested By:				
Testing Results				
Test	Result	Outcome	Cutoff	Notes
Chain of Custody Faults				
No Donor INITIALS on SEAL		*		
No Donor SIGNATURE on COC		*		
No COLLECTOR Signature on COC		*		
Screening Tests by IA				
URINE: Meth/Amphetamines		SEE BELOW	1000 ng/ml	
URINE: Benzodiazepines		negative	300 ng/ml	
URINE: Cocaine		negative	300 ng/ml	
URINE: Methadone		negative	300 ng/ml	
URINE: Opiates		negative	300 ng/ml	
URINE: THC		negative	50 ng/ml	
URINE: Oxycodone		negative	300 ng/ml	
URINE: Ethyl Glucuronide-ETG		negative	500 ng/ml	• ETG negative screen suggests alcohol was not consumed.
URINE: Buprenorphine/Suboxone		SEE BELOW	5 ng/ml	

- Faults noted on results- documenting errors in the collection process and chain of custody issues
- Statistical reporting on collection faults and quality steps- prior to testing

DETECTED VS POSITIVE- CONFIRMATION TESTING

Confirmation Tests				
Methadone by LC/MS/MS				
URINE: Methadone metabolite	1020 ng/ml	**POSITIVE	300 ng/ml	
URINE: Methadone	1375 ng/ml	**POSITIVE	300 ng/ml	
Opiates by LC/MS/MS				
URINE: Codeine	Not Detected	negative	300 ng/ml	
URINE: Morphine	Not Detected	negative	300 ng/ml	
URINE: Hydrocodone	380 ng/ml	**POSITIVE	300 ng/ml	
URINE: Hydromorphone	250 ng/ml	Detected	300 ng/ml	• Results consistent with Hydrocodone use. • Detected (drug fingerprint confirmed) above the laboratory limit of quantitation.

Drug Class	LOQ ng/mL	Cut Off ng/mL
THC	3.5	15
Cocaine	12	150
Amp/Meth	75	250/500
Opiates	25	300
Benzodiazepines	25	100
Tramadol	50	100
Buprenorphine	3.75	5

Ex. 25-299 ng/mL → “Detected”
≥300 ng/mL → “Positive”

LOQ = Lowest limit of quantitation

Below the established “cut off” but still absolutely present confirming the presence of drug in the sample

ORAL FLUID DRUG TEST FINAL REPORT



Flagstaff Lab
1760 E Route 66
Flagstaff, AZ 86004
855-895-8090

Denver Lab
1701 Chambers Rd, Unit J
Aurora, CO 80011
855-895-8090

Long Island Lab
789 Park Avenue
Huntington, NY 11743
855-895-8090

Tacoma Lab
2617 East L Street
Tacoma, WA 98421
855-895-8090

Worcester Lab
415 Main Street, 4th Floor
Worcester, MA 01608
855-895-8090

Specimen Information				
Donor Name:	DOE, JOHN	Collected:	10/01/2020 13:55:00	
Donor DOB:	12/25/2000	Received:	10/01/2020 10:55:11	
Accession:	0Z2000138	Reported:	10/01/2020 11:07:43	
COC:	AJ563535DC	Donor Other ID:		
Type (Matrix):	Oral Fluid	Donor Case:	5465321	
Client Code:	TRM			
Client:	TRM, INC. NEW HIRE			
Requested By:	JANE DOE			
Testing Results				
Test	Result	Outcome	Cutoff	Notes
Screening Tests by IA				
ORAL FLUID: Ethanol		negative	0.02 %	
ORAL FLUID: Meth/Amphetamines(50 ng/ml)		**POSITIVE	50 ng/ml	
ORAL FLUID: Benzodiazepine		negative	20 ng/ml	
ORAL FLUID: Methadone		negative	50 ng/ml	
ORAL FLUID: Cannabinoids/THC		negative	4 ng/ml	
ORAL FLUID: Cocaine		negative	20 ng/ml	
ORAL FLUID: Opiates		negative	40 ng/ml	
ORAL FLUID: Buprenorphine		negative	5 ng/mL	
ORAL FLUID: PCP		**POSITIVE	10 ng/ml	
Confirmation Tests				
Meth/Amphetamine Oral Fluid by LC/MS/MS				
Amphetamine-OF	48 ng/ml	**POSITIVE	10 ng/ml	
Methamphetamine-OF	136 ng/ml	**POSITIVE	10 ng/ml	
MDA-OF		negative	10 ng/ml	
MDMA-OF		negative	10 ng/ml	
MDEA-OF		negative	10 ng/ml	
PCP Oral Fluid by LC/MS/MS				
PCP-OF	31.2 ng/ml	**POSITIVE	4 ng/ml	
Additional Comments				
• Testing performed at Cordant Internal Accessions, , , . • Confirmation of a positive screen is recommended if legal action is anticipated. • Specimen received sealed and intact unless otherwise noted. • Report Released By: TJG - B.S., Certifying Scientist.				

- Mostly parent drugs, not metabolites
- Lower cut offs than urine
- Very recent use (1-3 days)

RESOURCES



TRAINING OPPORTUNITIES

- Visit our website for recorded trainings and webinars
 - Sign up for alerts to future live webinars available to the public
- Visit us on Youtube for quick 2–3-minute short trainings on most frequently asked questions

The screenshot shows the Cordant Health Solutions website with a navigation bar including links for INDUSTRIES, SOLUTIONS, TESTING OPTIONS, ABOUT US, RESOURCES, NEWS, and a CONTACT US button. Below the navigation bar, there are three featured training modules, each with a video thumbnail and a description.

Training Title	Instructor(s)	Description	Action
UNDERSTANDING ORAL FLUID TESTING AND RESULT INTERPRETATION	Cynthia Whiteman, MS, D-ABFT-FT	Saliva Testing 101 Learn more about saliva drug testing—how it is used, how to interpret test results, and why it may be a good fit for your practice or agency now and after the current health crisis.	WATCH »
ENSURING PATIENT ACCESS TO BUPRENORPHINE DURING THE COVID-19 PANDEMIC	Demetra Ashley, Pharmacy Compliance Officer Bob Mann, Senior VP & GM of Integrated Services	Ensuring Patient Access to Buprenorphine During the COVID-19 Pandemic Cordant's pharmacy solution can help ensure patients have access to buprenorphine during the COVID-19 pandemic.	WATCH »
MENTAL HEALTH MEDICATION ADHERENCE - WHAT ARE YOU MISSING?	Cynthia Whiteman, MS, D-ABFT-FT	Mental Health Adherence Testing - What Are You Missing? Whether you work in criminal justice or drug courts, addiction treatment or pain management, if you're not drug testing for mental health medication adherence, what are you missing?	WATCH »

<https://cordantsolutions.com/cordant-videos/>

"I love all of the Webinars, this is the only trusted place I can get training for my job. Please do more!"

"These are great and very informational."

"All information related to drug screen collection/interpretation is useful and the more recent the better!"

EDUCATIONAL RESOURCES

| ☎ 844-848-5955 | ✉ info@cordanths.com

Login

Cordant

Health Solutions®

[INDUSTRIES](#)
[SOLUTIONS](#)
[TESTING OPTIONS](#)
[ABOUT US](#)
[RESOURCES](#)
[NEWS](#)
[CONTACT US](#)

Barbiturates

Barbiturates refers to a class of drugs that affect the central nervous system as a depressant. Common barbiturates include phenobarbital and butabarbital. Fill out the form now to download the pdf.

[READ MORE »](#)

Benzodiazepine Metabolism

Benzodiazepines belong to the sedative-hypnotic group of drugs that possess selective central nervous system, CNS depressant activity.

[READ MORE »](#)

Benzodiazepines

Benzodiazepines refers to a class of drugs that affects the central nervous system as a depressant. Common benzodiazepines include lorazepam, clonazepam and diazepam.

[READ MORE »](#)

<https://cordantsolutions.com/drug-education-resource-library/>


DRUG DETECTION TIMES

Drug Class and Drug	Common Prescription and Brand Name(s)	Drug(s) Detected	Predicted Maximum Detection Times (Days)	
Alcohol/Alcohol Biomarkers			Urine	Oral Fluid
Alcohol	N/A	Ethyl Alcohol	<1	<1
Alcohol Biomarkers	N/A	Ethyl Glucuronide/Ethyl Sulfate	1-3	<1
Amphetamines/Stimulants			Urine	Oral Fluid
Amphetamine/ Dextroamphetamine	Adderall, Benzedrine, Dexedrine	Amphetamine	1-5	1-2
Methamphetamine	Desoxyn, Metabolite of Didrex, Select Generic Medicated Vapor Inhalers	Amphetamine*, Methamphetamine	1-5	1-3
Methylphenidate	Concerta, Daytrana, Focalin, Metadate, Methylin, Ritalin	Methylphenidate, Ritalinic Acid*	1-2	1
Antianxiety/Benzodiazepines			Urine	Oral Fluid
Alprazolam	Xanax	Alpha-Hydroxyalprazolam*, Alprazolam	1-4	1-2
Chlordiazepoxide	Librium	Nordiazepam*, Oxazepam*	1-10	1-2
Clonazepam	Klonopin	7-Aminoclonazepam*, Clonazepam	1-4	1-2
Diazepam	Valium	Diazepam, Nordiazepam*, Oxazepam*, Temazepam*	1-7	1-2
Lorazepam	Ativan	Lorazepam	1-7	1-2
Oxazepam	Serax	Oxazepam	1-7	1-2
Temazepam	Restoril	Oxazepam*, Temazepam	1-4	1-2
Anticonvulsants/Antiepileptics			Urine	Oral Fluid
Gabapentin	Gralise, Neurontin	Gabapentin	1-4	1-2
Pregabalin	Lyrica	Pregabalin	1-4	1-2

Drug Detection Chart provides helpful cross reference for testing information

- Drug and Drug Class
- Common Brand and Prescription Names
- Drugs/Metabolites Detected In Our Testing
- Average Maximum Detection Window for each drug in Urine and/or Oral fluid

BENZODIAZEPINE METABOLISM




Benzodiazepines belong to the sedative-hypnotic group of drugs that possess selective CNS depressant activity. Sedative properties refer to calming or tranquilizing effects of the drug. Hypnotic properties are the ability to induce sleep. The sedative - hypnotic properties of the various benzodiazepines overlap to a degree and are dose-dependent. The sedative properties are used to treat anxiety disorders, seizures, alcohol withdrawal, as muscle relaxants and as preanesthetic agents. The hypnotic and anesthetic properties require much higher doses compared to sedative effects. The indicated uses on page 2 are general uses and may vary depending on individual clinical situations.

Detection times for benzodiazepine are approximate and vary depending on dosage, frequency of dosing and the length of time that a patient has been using a particular drug. Oxazepam may be detected for several weeks after administration of long-acting drugs such as diazepam or chlordiazepoxide or much shorter detection times after administration of intermediate-acting drugs like temazepam or oxazepam.

OPIOID METABOLISM

This reference chart details the metabolism of several key opioids to aid in the interpretation of toxicology testing results. An arrow pointing towards a specific drug/metabolite indicates a potential source, while an arrow pointing away from a specific drug/metabolite indicates a metabolic product.


For example, pharmaceutical morphine, codeine, heroin and poppy seeds are all possible sources of a morphine positive test result. Approximately 15-20% of codeine is metabolized into morphine. Heroin itself has an extremely short half life; it is metabolized into 6-acetylmorphine (6-AM). 6-AM is a very unique metabolite of heroin; thus if a specimen tests positive for 6-AM, most likely it is due to heroin use. However, 6-AM is also further metabolized into morphine. Poppy seeds also contain morphine and a very small amount of codeine. Reference the table on the reverse side for more information.



SUPPORT SERVICES



**CUSTOMER
SERVICE**



**TOXICOLOGY
INTERPRETATION**



**LITIGATION
SUPPORT**



**REFERRAL
SUPPORT**

**Dedicated Indiana DCS
Support Hotline**

866-217-0862

ORDERING COLLECTION SUPPLIES:

Select Option 3 for the Client Services Team or e-mail Cordantclientservices@cordanth.com to order supplies. Reference your account number that will begin with CS.